CE 0459*

Universal Dental Adhesive

Filled Universal Dental Adhesive

Instructions for Use

U.S. Patents: 5,749,733 & 5,658,963

IN-130R7

Rev. 10/18

ONE-STEP* PLUS Filled Universal Dental Adhesive
GENERAL INFORMATION

ONE-STEPONE-STEP PLUS (OSOs®) are technological breakthroughs from BISCO. OSOs® are 5th generation (2-step) Universal Adhesives. **ONE-STEP PLUS** is 8.5% filled with a reflowing proprietary filler with an average particle size of 1µm. **OSOs®** are designed for most composites to dentin, enamel, cast metals, treatment porcelain and set amalgam. Because of its unique chemistry, and its minimal and uniform film thickness, this system is also designed for microleakage when used with a conventional medium as well as amalgam bonding. **OSOs®** are truly universal adhesives that have been designed to work with light-cured, self-cured and dual-cured composite materials.

Indications for use:

- All Direct Composite Restorations
- Direct Shrinkage Composite Restorations
- Desensitization of Crown Preparations Prior to Impression Making/Provisionalization
- Indirect Restorative Cementation, Metal, Porcelain and Composite
- Adhesive (Bonded) Amalgam Restorations
- Fresh Amalgam to Old Amalgam
- Composite Core Build-Ups
- Porcelain Repairs
- Composite to Composite (Repair)
- Composite to Metal/Steel Amalgam (Direct Veneering)
- Root Desensitization
- Endodontic Post Bonding

Warnings:

- Highly Flammable
- When using dental adhesives, saliva contamination will seriously compromise dentin bonding.
- The phosphoric acid in the etchant is a severe eye and skin irritant. Hryngy may irritate the etchant is allowed to remain on the skin or mucosa for extended periods of time. If accidentally splashed into the eyes, flush with copious amounts of water and seek medical attention immediately. In case of contact with other tissues, rinse immediately with plenty of water for several minutes.
- A bitter taste is essential to the success of a porcelain repair. Measure the composition of porcelain or metal is dentral.

Cautions:

- Cross-contamination: Product may contain fibers that are used for one time use. Dispose of used or contaminated unit-dose components and accessories. Do not clean, disinfect or reuse.

Precautions:

- For **ONE-STEP PLUS**, shake bottle well to redistribute any filler particles that may have settled. If needed, invert bottle and tap on hard surface to help agglutator.
- For optimal performance, use when **OSOs®** are at room temperature.
- Avoid contact with the skin, unpolymerized (meth)acrylate resins may cause skin sensitization in susceptible persons. In case of contact, wash skin with soap and water.
- Refer to HELPFUL HINTS section for additional technique-related precautions.
- Use individual component labels for specific expiration dates.
- Safety data sheet available on request.
- Safety data sheet available on request.

HELPFUL HINTS

OSOs® can be used with bases or liners of metallic, resin-modified glass ionomer or glass ionomer.

Precautions:

- ZOC bases may inhibit the set of the adhesive or composite.
- Covering the dentin with a base or liner prior to adhesive application may compromise dentin bonding.
- For optimal performance, use when **OSOs®** are at room temperature.
- Avoid contact with the skin, unpolymerized (meth)acrylate resins may cause skin sensitization in susceptible persons. In case of contact, wash skin with soap and water.
- Refer to HELPFUL HINTS section for additional technique-related precautions.
- Use individual component labels for specific expiration dates.
- Safety data sheet available on request.
- Safety data sheet available on request.

3 OSOs® adhesives contain a monomer, BPDM, dissolved in an acetone solvent. Once 2 coats have been applied to the dentin surface, though not being air dried for up to 10 seconds should be done to remove excess solvent. Begin with a gentle stream initially and increase to a fairly strong stream of air.

• Inadequate air drying will not completely evaporate the solvent and may lead to reduced bond strength and post-operative sensitivity.

Proper air drying is extremely important with all bonding procedures but particular attention should be given to air drying of inlay or onlay preparations. Care should be given to avoid pooling of adhesive in the internal line angles. A stream of air should be directed at all line angles to properly thin adhesive.

• Excessive pooling could interfere with the setting of the definitive restoration.

To increase adhesion to exposed metal, a sealant is highly recommended for surface roughening. Diamond bur abrasion may not provide enough surface area for proper bonding.

6. Always lubricate the matrix band prior to amalgam bonding procedures to prevent adhesion of amalgam to the band. Hard inlay wax or a crown waxes will. Copalite or vaseline are not always available and therefore, not recommended.

7. Adhesive (Bonded) Amalgam Restorations require the use of RESINOMER™, a dual-cured Amalgam Bonding/Luting System (sold separately). RESINOMER bonds tenaciously to **OSOs®** as well as amalgam. Mechanical adhesion of the amalgam is condensed into the uncured RESINOMER, while chemical adhesion occurs due to the diarylmethoxy dimethacrylate (DSDM).

INSTRUCTIONS FOR USE

- Isolate tooth and prepare the cavity in a conservative manner. Clean the entire surface with a slurry of pumice and a cleaner, such as CAVITY CLEANSER™, or water.
- Etch enamel and dentin using an etchant, such as UNI-ETCH™ wBAC or ETCH-37™™ wBAC for 15 seconds.
- Prepared enamel: Rinse thoroughly. Remove excess water by blotting the surface with a Foam Pellet, leaving the dentin and enamel visibly moist with a shiny surface.
- Inadequate air-drying will not completely evaporate the solvent and may lead to reduced bond strength and post-operative sensitivity.
- Proper air-drying is extremely important with all bonding procedures but particular attention should be given to air drying of inlay or onlay preparations. Care should be given to avoid pooling of adhesive in the internal line angles. A stream of air should be directed at all line angles to properly thin adhesive.
- Excessive pooling could interfere with the setting of the definitive restoration.
- To increase adhesion to exposed metal, a sealant is highly recommended for surface roughening. Diamond bur abrasion may not provide enough surface area for proper bonding.
- Always lubricate the matrix band prior to amalgam bonding procedures to prevent adhesion of amalgam to the band. Hard inlay wax or a crown waxes will. Copalite or vaseline are not always available and therefore, not recommended.

7. Adhesive (Bonded) Amalgam Restorations require the use of RESINOMER™, a dual-cured Amalgam Bonding/Luting System (sold separately). RESINOMER bonds tenaciously to **OSOs®** as well as amalgam. Mechanical adhesion of the amalgam is condensed into the uncured RESINOMER, while chemical adhesion occurs due to the diarylmethoxy dimethacrylate (DSDM).

INSTRUCTIONS FOR USE

- Isolate tooth and prepare the cavity in a conservative manner. Clean the entire surface with a slurry of pumice and a cleaner, such as CAVITY CLEANSER™, or water.
- Etch enamel and dentin using an etchant, such as UNI-ETCH™ wBAC or ETCH-37™™ wBAC for 15 seconds.
- Prepared enamel: Rinse thoroughly. Remove excess water by blotting the surface with a Foam Pellet, leaving the dentin and enamel visibly moist with a shiny surface.
- Inadequate air-drying will not completely evaporate the solvent and may lead to reduced bond strength and post-operative sensitivity.
- Proper air-drying is extremely important with all bonding procedures but particular attention should be given to air drying of inlay or onlay preparations. Care should be given to avoid pooling of adhesive in the internal line angles. A stream of air should be directed at all line angles to properly thin adhesive.
- Excessive pooling could interfere with the setting of the definitive restoration.
- To increase adhesion to exposed metal, a sealant is highly recommended for surface roughening. Diamond bur abrasion may not provide enough surface area for proper bonding.
- Always lubricate the matrix band prior to amalgam bonding procedures to prevent adhesion of amalgam to the band. Hard inlay wax or a crown waxes will. Copalite or vaseline are not always available and therefore, not recommended.

NOTE: Air dry the surface thoroughly to ensure complete evaporation of the acetone solvent. Once 2 coats have been applied to the dentin surface, though not being air dried for up to 20-25 seconds. Blot the preparation dry with a Foam Pellet or gently air dry for 1-2 seconds to remove excess moisture.

4. For bottle delivery:
Precautions:

- When using **ONE-STEP PLUS**, shake bottle well to ensure the mixing element inside the bottle is audible. If needed, invert bottle and tap on hard surface to help agglutator.
- For optimal performance, use when **OSOs®** are at room temperature.
- Avoid contact with the skin, unpolymerized (meth)acrylate resins may cause skin sensitization in susceptible persons. In case of contact, wash skin with soap and water.
- Refer to HELPFUL HINTS section for additional technique-related precautions.
- Use individual component labels for specific expiration dates.
- Safety data sheet available on request.
- Safety data sheet available on request.

OPTIONAL FOR CLASS V RESTORATIONS: Apply 1-2 additional coats of **OSOs®** to the entire internal surface of the preparation, agglating slightly for 10 seconds. Thoroughly air dry and light-cure for 10 seconds.

7. Proceed with restorative treatment.

Research indicates a thicker layer of adhesive will help dissipate the effects of tooth flexure while providing a wai air-inhibited resin layer to help adapt dry, viscous composites.**

WARNING: Store at room temperature (20°C/68° F - 25°C/77° F).

STORAGE: BISCO, Inc. recognizes its responsibility to reproduce products it proven to be defective. BISCO, Inc. does not accept liability for any damages or loss, either direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection with its use.

1. Kancz, J. One-Step bond to enamel and dentin. Am. J Dent 1997;10:1 5-8.
2. Kemp-Scholtz CM, Davidson CL. Marginal integrity related to bond strength and strain capacity of composite resin restorative systems. J Prosthet Dent 1995;64:684-694.
3. Van Meerbeek, B. Dentine Adhesion: Morphological, physico-chemical and Clinical Aspects. Thesis, Leven (Belgium) 1993.

ONE-STEP and UNI-ETCH are registered trademarks of BISCO, Inc. RESINOMER, CAVITY CLEANSER, ETCH-37 and AQUA-PRAP are trademarks of BISCO, Inc. Microbrush is a trademark of Microbrush, Int'l.

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