CORE-FLO™ DC Lite
Low Viscosity Core Build-Up

Dual-Cured

BISCO, Inc.
1100 W. Irving Park Road
Schamburg, IL 60193
U.S.A.
1-800-247-3368
1-847-534-6000
www.bisco.com

BISICO France
208, allée de la Coudoulette
13680 Lançon de Provence
France
33-4-90-42-92-92

IN-218R7
Rev.10/20

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed Dentist.

GENERAL INFORMATION

CORE-FLO DC Lite is a flowable dual-cured, radiopaque core build-up composite, specifically designed to cement posts and build-up the core structure in one easy application. **CORE-FLO DC Lite** is also an ideal dentin replacement material. **CORE-FLO DC Lite** comes in an auto mix dual-syringe and is available in both Natural/A1 & Opaque White shades. **CORE-FLO DC Lite** is intended for use with a bonding agent that is compatible with dual-cured materials, such as BISCO adhesives.

- Indications for use:**
The principle indications for use of **CORE-FLO DC Lite** are:
1. Core build-up material
2. Cementation of posts or pins
3. Dentin replacement material
Intended user: licensed dentist.

- Warnings:**
- Contamination by saliva, blood, or water should be avoided during these procedures. When using dental adhesives, contamination will compromise dentin bonding and may result in decreased longevity of the restoration.
 - A rubber dam is recommended for all indications where contamination may be possible.

- Cautions:**
- Cross-contamination: Product may contain items that are designed to be used with dental syringes or contaminated accessories. Do not clean, disinfect or reuse.
 - Use of commonly available hygienic protective covering to avoid any contamination of syringes during treatment is recommended. If syringe becomes contaminated, discard. Do not clean or disinfect.

- Precautions:**
- Avoid contact with the skin; unpolymerized (meth) acrylate resins may cause skin sensitization in susceptible persons. In case of contact, wash skin with soap and water.
 - To prevent polymerization or clogging, the dual-syringe mixing tip should be left in place until the next application. If used intraorally, wipe mixing tip with disinfectant prior to storage.
 - Unlike Bisco adhesives, not all one bottle total-etch or self-etch systems bond to self-cured or dual-cured materials. Refer to the manufacturer's instructions for any bonding with these materials.
 - If refrigerated, product should be brought to room temperature before use.
 - See individual component labels for specific expiration dates.
 - Safety data sheet available on request.
 - Safety data sheet available at www.bisco.com.
 - Instructions for Use available upon request.
 - Instructions for Use available at www.bisco.com.

TECHNICAL INFORMATION
Air room temperature (20°C/68°F - 25°C/77°F)
Working Time: Minimum 1 minute
Setting Time: Maximum 7 minutes
Working time and Setting time may vary based on storage conditions, temperature, humidity, etc.

INSTRUCTIONS FOR USE

CORE MATERIAL MIXING PROCEDURES:
1. Remove the syringe cap or mixing tip from the dual-syringe.
2. Dispense a small amount of material onto a mixing pad to eliminate any voids in each chamber of the dual-syringe.
3. Attach a mixing tip to the dual-syringe by lining up the key and keyway. Then turn the brown locking mechanism clockwise.
4. Pressing the plunger will mix and dispense the **CORE-FLO DC Lite**. (See "TECHNICAL INFORMATION" Section for Working and Setting Time information.)

1. CORE BUILD-UP USING CORE-FLO DC Lite
(1) Isolate the tooth and prepare the tooth according to the post manufacturer's instructions; we recommend the use of any BISCO adhesive.
(2) If using a matrix band: place a matrix band around the prepared tooth. Place the end of the intra-oral tip directly on to the preparation and dispense **CORE-FLO DC Lite** into the preparation. Push the intra-oral tip into the preparation and dispense **CORE-FLO DC Lite**. Place the loaded core onto the prepared tooth and hold in place, light cure for 40 seconds.
(4) Prepare and finish the core in the usual manner.

2. POST CEMENTATION AND CORE BUILD-UP TECHNIQUE USING CORE-FLO DC Lite
(1) Isolate tooth and prepare the tooth according to the post manufacturer's instructions.
(2) Apply an adhesive according to the manufacturer's instructions; we recommend the use of any BISCO adhesive. Prior to light curing, blot the canal dry with paper points until the paper point returns dry from the canal and gently air dry the adhesive to evaporate the solvent.
NOTE: The paper point step is important for removal of any pooled adhesive so that there is no interference with the setting of the post.
(3) Coat the surface of the post with an adhesive according to the manufacturer's instructions. This will increase the wettability of the post.
(4) Using the root canal mixing tip inject **CORE-FLO DC Lite** into the canal. Start at the apical end of the canal and withdraw the tip slightly while keeping the tip in the cement to minimize voids.
(5) Coat the apical end of the post with **CORE-FLO DC Lite**. Seal the post into the canal gently and maintain firm pressure for 20-30 seconds until the post is seated.
(6) Continue to express **CORE-FLO DC Lite** around the post and build-up to the core to the desired height and shape.
(7) Remove the excess material and light cure for 40 seconds by placing the light tip as close as possible to the post and core build-up.
(8) Prepare and finish the core in the usual manner.

3. DENTIN REPLACEMENT MATERIAL
(1) Isolate tooth and prepare in a conservative manner. Rinse and dry.
(2) Apply an adhesive according to the manufacturer's instructions; we recommend the use of any BISCO adhesive.
(3) Express **CORE-FLO DC Lite** only to the level of the dentin-enamel junction.
(4) Light cure the **CORE-FLO DC Lite** for 10 seconds. Continue to firmly pack a light-cured composite, building incrementally until cavity preparation is filled. Light cure according to manufacturer's instructions.

HYGIENE: To prevent polymerization or clogging, the dual-syringe mixing tip should be left in place until the next application. Use of commonly available hygienic protective covering to avoid any contamination of syringes during treatment is recommended. If used intraorally, wipe mixing tip with disinfectant prior to storage.
DISPOSAL: Refer to community provisions relating to waste. In their absence, refer to national or regional provisions relating to waste.
STORAGE: Store kit at room temperature (20°C/68°F - 25°C/77°F). See individual component labels for specific expiration dates.
WARRANTY: BISCO, Inc. recognizes its responsibility to replace products if proven to be defective. BISCO, Inc. does not accept liability for any damages or loss, either direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

*CORE-FLO is a trademark of BISCO, Inc.
To Reorder Accessories:
X-81252P – Dual-Syringe Static Mixing Tips (30)
X-81257P – Dual-Syringe Mixing and Intraoral Tips (30)
X-81267P – Dual-Syringe Mixing and Root Canal Tips (30)

NOTICE: Report any serious incident to BISCO and the competent authority of the EU member state in which the user and/or patient is established.

3. DENTIN REPLACEMENT MATERIAL
(1) Isolate tooth and prepare in a conservative manner. Rinse and dry.
(2) Apply an adhesive according to the manufacturer's instructions; we recommend the use of any BISCO adhesive.
(3) Express **CORE-FLO DC Lite** only to the level of the dentin-enamel junction.
(4) Light cure the **CORE-FLO DC Lite** for 10 seconds. Continue to firmly pack a light-cured composite, building incrementally until cavity preparation is filled. Light cure according to manufacturer's instructions.

HYGIENE: To prevent polymerization or clogging, the dual-syringe mixing tip should be left in place until the next application. Use of commonly available hygienic protective covering to avoid any contamination of syringes during treatment is recommended. If used intraorally, wipe mixing tip with disinfectant prior to storage.
DISPOSAL: Refer to community provisions relating to waste. In their absence, refer to national or regional provisions relating to waste.
STORAGE: Store kit at room temperature (20°C/68°F - 25°C/77°F). See individual component labels for specific expiration dates.
WARRANTY: BISCO, Inc. recognizes its responsibility to replace products if proven to be defective. BISCO, Inc. does not accept liability for any damages or loss, either direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

*CORE-FLO is a trademark of BISCO, Inc.
To Reorder Accessories:
X-81252P – Dual-Syringe Static Mixing Tips (30)
X-81257P – Dual-Syringe Mixing and Intraoral Tips (30)
X-81267P – Dual-Syringe Mixing and Root Canal Tips (30)

NOTICE: Report any serious incident to BISCO and the competent authority of the EU member state in which the user and/or patient is established.

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AVFOER: Rådplege de varemærke tilhørende BISCO, Inc.