BISCO ETCHANTS


Semi-Gel

Phosphoric Acid Etchant

UNI-ETCH® (32%) W/BAC

ETCH-37™ (37%) W/BAC

Instructions for Use

BISCO, Inc.
1100 W. Irving Park Road
Schaumburg, IL 60193
U.S.A.
1-800-247-3368
1-847-534-6000
www.bisco.com

ISO 9001

BISCO France
208, allée de la Coudoulette
13680 Lançon de Provence
France
33-4-90-42-92-92

ISO 9001

BISOC France
208, allée de la Coudoulette
13680 Lançon de Provence
France
33-4-90-42-92-92

IN-03919

Rev. 5/20

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed Dentist.

| | |
|-----------|---|
| en | ETCHANTS UNI-ETCH® (32%) UNI-ETCH™ (32%) W/BAC ETCH-37™ (37%) ETCH-37™ (37%) W/BAC |
| | Phosphoric Acid Etchant |
| | GENERAL INFORMATION |

BISCO ETCHANTS are phosphoric acid (H₃PO₄) intended for use on enamel, dentin and as a descaling agent on restorative substrates. BISCO ETCHANTS will effectively remove the smear layer, etch and demineralize enamel and dentin, producing the necessary microretentive surface for successful bonding. The etchant rinses away cleanly and quickly, leaving no residue. BISCO ETCHANTS are also available with benzalkonium chloride (BAC), an antimicrobial agent. Research confirms the development of residual bacterial inhibition zones of up to 7 mm, when challenged (in-vitro) with Actinomyces viscosus. **NOTE:** The inclusion of benzalkonium chloride has not been shown to correlate with a reduction in secondary decay in patients. In-vitro clinical studies to evaluate the effects of benzalkonium chloride on oral bacteria are still not being performed.

Intended user: licensed dentist.

Indications For Use:

- Etching dentin and enamel
- Selective enamel etching
- Cleaning agent on dental restorative materials

Warnings:

- The phosphoric acid in the etchant is a severe eye and skin irritant. Injury may result if the etchant is allowed to remain on the skin or mucosa for extended periods of time. In case of contact with other tissues, rinse immediately with plenty of water for several minutes.
- If accidentally splashed into the eye, flush with copious amounts of water and seek medical attention immediately.

Precautions:

Cross-contamination:

- Product may contain items that are designed for one time use. Dispose of used or contaminated tips. Do not clean, disinfect or reuse.
- Use of commonly available hygienic protective covering to avoid any contamination of syringes during treatment is recommended. If syringe becomes contaminated, discard. Do not clean or reuse.

Keep away from direct sunlight. Prolonged exposure to sunlight could cause a discoloration of the etchant; this will not compromise the efficacy of the product.

After the application of Etchant on a mixing pad or a glass slab, BEFORE use on patients so as to become familiar with the pressure required to dispense the Etchant from the tip. Never use introrally before checking flow on tip. Never use force on plungers. Start from introrally. This could result in a substantial uncontrolled discharge of material. If resistance to flow is experienced, DO NOT proceed. Replace with new tip and verify flow before use on patient.

Prophy pastes containing oil or fluoride should not be used as these additives may interfere with etching.

- See individual component labels for specific expiration dates.

A safety data sheet available on request.

For more information visit www.bisco.com.

Instructions for Use available upon request.

INSTRUCTIONS FOR USE

- Etching dentin and enamel**
 - Remove white cap from etchant syringe.
 - Attach the light blue 25 gauge tip.
 - Verify flow of etch before intraral application. Dispense a small portion of gel onto a pad or gauze to familiarize yourself with its viscosity.
 - If the surface to be etched is not freshly cut, clean with a pump-water syringe. Prophy pastes containing oil or fluoride should not be used as these additives may interfere with etching. Rinse well.
 - Isolate the preparation.
 - Place etchant directly onto the dentin and enamel surfaces to be etched.
 - Use wettable tip piece for 15 seconds.
 - Thoroughly rinse with water for 5 seconds with suction.
 - Immediately proceed with the bonding procedure.
 - Discard all single use components.
- Selective enamel etching**
 - Remove white cap from etchant syringe.
 - Attach the light blue 25 gauge tip.
 - Verify flow of etch before intraral application. Dispense a small portion of gel onto a pad or gauze to familiarize yourself with its viscosity.
 - If the surface to be etched is not freshly cut, clean with a pump-water syringe. Prophy pastes containing oil or fluoride should not be used as these additives may interfere with etching. Rinse well.
 - Isolate the preparation.
 - Place etchant directly onto the surfaces to be etched. Leave etchant in place for 10 seconds.
 - Thoroughly rinse with water for 5 seconds with suction.1. selectively etching enamel, suction prior to rinsing to avoid etchant flowing on to dentin.
 - Disinfect area around preparation with the bonding procedure.
 - Discard all single use components.
- Cleaning agent on dental restorative materials**
 - BISCO ETCHANTS may be used on restorative substrates to effectively remove organic contaminants. Apply etchant to the surface for 30 seconds and rinse thoroughly with water. **NOTE:** Unprimed zirconia surfaces should NOT be cleaned with phosphoric acid. Use a restoration cleaner such as ZiClean to remove contaminants.

HYGIENE: Re-cap syringes with the lock cap.

DISPOSAL: Refer to community provisions relating to waste. In their absence, refer to national or regional provisions relating to waste.

STORAGE: Store at room temperature (20°C/68°F - 25°C/77°F). Store away from direct sunlight.

WARRANTY: BISCO, Inc. recognizes its responsibility to replace products if proven to be defective. BISCO, Inc. does not accept liability for any damages or loss, either direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

- UNI-ETCH is a registered trademark of BISCO, Inc.
- PORCELAN PRIMER is manufactured by BISCO, Inc.
- ETCH-37, BIS-SILANE and Z-PRIME are trademarks of BISCO, Inc.

NOTICE: Report any serious incident to BISCO and the competent authority of the state in which the user or patient is established.

TO REORDER:

X-80611N – 125 Light Blue Disposable Tips (25 Gauge)
X-80608N – 50 Dark Blue Disposable Tips (22 Gauge)

- Kanca, J. J. "Etchant Composition and Bond Strength to Dentin". Am J Dent 1993;6:162-164.

- Dr. Daniel Chan, University of Texas Health Science Center at San Antonio Dental School,Residual Effect of 1 and 2% Benzalkonium Chloride Incorporated into an Etchant on the Susceptibility of Actinomyces viscosus T14V 1993.

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WARRANTY: We acknowledge our responsibility to replace products if proven defective. BISCO, Inc. does not assume liability for damages or loss, either direct or consequential, resulting from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

- Kanca, J. J. "Etchant Composition and Bond Strength to Dentin". Am J Dent 1993;6:162-164.

- Dr. Daniel Chan, University of Texas Health Science Center at San Antonio Dental School,Residual Effect of 1 and 2% Benzalkonium Chloride Incorporated into an Etchant on the Susceptibility of Actinomyces viscosus T14V 1993.

| | |
|-----------|--|
| de | ÄTZMITTEL UNI-ETCH® (32%) UNI-ETCH™ (32%) W/BAC ETCH-37™ (37%) ETCH-37™ (37%) W/BAC |
| | Phosphorsäure-Ätzmittel |
| | ALLGEMEINE INFORMATIONEN |

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- Remove white cap from the etchant syringe.
- Attach the light blue 25-gauge tip.
- Verify flow of etch before intraral application. Dispense a small portion of gel onto a pad or gauze to familiarize yourself with its viscosity.
- If the surface to be etched is not freshly cut, clean with a pump-water syringe. Prophy pastes containing oil or fluoride should not be used as these additives may interfere with etching. Rinse well.
- Isolate the preparation.
- Place etchant directly onto the surfaces to be etched. Leave etchant in place for 15 seconds.
- Thoroughly rinse with water for 5 seconds with suction.
- Immediately proceed with the bonding procedure.
- Discard all single use components.

Cleaning agent on dental restorative materials

BISCO ETCHANTS may be used on restorative substrates to effectively remove organic contaminants. Apply etchant to the surface for 30 seconds and rinse thoroughly with water.

NOTE: Unprimed zirconia surfaces should NOT be cleaned with phosphoric acid. Use a restoration cleaner such as ZiClean to remove contaminants.

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