BISCO's ALL-BOND UNIVERSAL is a light-cured single-component dental bonding agent that combines etching, priming and bonding in one bottle. ALL-BOND UNIVERSAL is an ethanol/water-based dental adhesive which bonds to dentin and cut and un-cut enamel. As a universal adhesive, ALL-BOND UNIVERSAL has been designed to be fully compatible with light-cured, self-cured and dual-cured composite materials. ALL-BOND UNIVERSAL can be used for the bonding of both direct and indirect restorations. ALL-BOND UNIVERSAL may be used with or without phosphoric acid etchant.

- Indications for Use: ALL-BOND UNIVERSAL Adhesive is used for:
1. ALL Direct Restorations (resin-based composite, resin-modified glass ionomers, core build-ups, etc.)
 2. ALL Indirect Restorations (metal, glass ceramics, zirconia/alumina, etc.)
 3. Desensitization/Sealing of Tooth
 4. Intraoral Repair (chipped porcelain, additions to direct restorations, etc.)
 5. Protective Varnish for Glass Ionomer Based Fillings

- Contraindications:
- Patients who have a history of severe allergic reactions to methacrylate resins.

- Warnings:
- Avoid splashing into the eyes. If ALL-BOND UNIVERSAL comes into contact with the eyes, flush with copious amounts of water and seek medical attention.
 - The phosphoric acid in the etchant is a severe eye and skin irritant. Injury may result if the etchant is allowed to remain on the skin or mucosa for extended periods of time. If accidentally splashed into the eye, flush with copious amounts of water and seek medical attention immediately. In case of contact with other tissues, rinse immediately with plenty of water for several minutes.
 - When using dental adhesives, saliva contamination will compromise dentin bonding and may result in decreased longevity of the restoration. A rubber dam is recommended for all direct restorations and other indications where contamination may be possible.
 - Highly Flammable.

- Cautions:
- Cross-contamination: Product may contain items that are designed for one time use. Dispose of used or contaminated accessories. Do not clean, disinfect or reuse.

- Precautions:
- ALL-BOND UNIVERSAL is a light-curable material. Use ALL-BOND UNIVERSAL immediately once it is dispensed in the mixing well. Prolonged exposure to air and light may lead to solvent evaporation and thickening of the adhesive.
 - To avoid evaporation of ALL-BOND UNIVERSAL, keep container tightly closed.
 - Avoid contact with the skin; unpolymerized (meth)acrylate resins may cause skin sensitization in susceptible persons. In case of contact, wash skin with soap and water.
 - If ZOE containing materials have been used on the tooth structure within 2-3 weeks, adequately clean the tooth surface with an abrasive prior to adhesive use.
 - The emitting tip of the dental curing light should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light-cured, it is advisable to divide the area into several sections and light cure each section separately.
 - Low curing light intensity causes poor adhesion.
 - See individual component labels for specific expiration dates.
 - Safety data sheet available on request.
 - Safety data sheet available at www.bisco.com.

INSTRUCTIONS FOR USE

- 1. Direct Restorations**

SELF-ETCHING TECHNIQUE: Prepare cavity. Wash thoroughly with water spray. Use an absorbent pellet or high volume evacuation (HVE) for 1-2 seconds, to remove excess water (do not desiccate). **NOTE:** BISCO's SELECT HV™ ETCH may be used on cut or uncut enamel to increase mechanical retention. This etchant was designed for precise selective placement. **TOTAL-ETCHING TECHNIQUE:** Prepare cavity. Wash thoroughly with water spray. Etch enamel and dentin using an etchant for 15 seconds. Rinse thoroughly. Remove excess water by blotting the surface with an absorbent pellet or high volume evacuation (HVE) for 1-2 seconds, leaving the preparation visibly moist.

 1. Dispense 1-2 drops of ALL-BOND UNIVERSAL into a clean well. Replace the cap onto the bottle promptly to avoid evaporation of solvent.
 2. Apply two separate coats of ALL-BOND UNIVERSAL, scrubbing the preparation with a microbrush for 10-15 seconds per coat. Do not light cure between coats.
 3. Evaporate excess solvent by thoroughly air-drying with an air syringe for at least 10 seconds; there should be no visible movement of the adhesive. The surface should have a uniform glossy appearance; otherwise, apply an additional coat of ALL-BOND UNIVERSAL and repeat Steps 2 and 3.
 4. Light cure for 10 seconds.
 5. Continue with placement of the restorative material according to the manufacturer's instructions.
- 2. Indirect Restorations**

Tooth Preparation: **SELF-ETCHING TECHNIQUE:** Prepare cavity. Wash thoroughly with water spray. Use an absorbent pellet or high volume evacuation (HVE) for 1-2 seconds, to remove excess water (do not desiccate). **NOTE:** BISCO's SELECT HV ETCH may be used on cut or uncut enamel to increase mechanical retention. This etchant was designed for precise selective placement. **TOTAL-ETCHING TECHNIQUE:** Prepare cavity. Wash thoroughly with water spray. Etch enamel and dentin using an etchant for 15 seconds. Rinse thoroughly. Remove excess water by blotting the surface with an absorbent pellet or high volume evacuation (HVE) for 1-2 seconds, leaving the preparation visibly moist.

 1. Dispense 1-2 drops of ALL-BOND UNIVERSAL into a clean well. Replace the cap onto the bottle promptly to avoid evaporation of solvent.
 2. Apply two separate coats of ALL-BOND UNIVERSAL, scrubbing the preparation with a microbrush for 10-15 seconds per coat. Do not light cure between coats.
 3. Evaporate excess solvent by thoroughly air-drying with an air syringe for at least 10 seconds; there should be no visible movement of the adhesive. The surface should have a uniform glossy appearance; otherwise, apply an additional coat of ALL-BOND UNIVERSAL and repeat Steps 2 and 3.
 4. Light cure for 10 seconds.
 5. Continue with placement of the restorative material according to the manufacturer's instructions.

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed Dentist. * B-7202P-AS is not CE marked.

BISCO, Inc. 1100 W. Irving Park Rd. Schaumburg, IL 60193 U.S.A. 1-847-534-6000 1-800-247-3368

EC REP BISCO France 208, allée de la Coudoulette 13680 Lançon De Provence France 33-4-90-42-92-92

BISCO, INC. 1100 W. Irving Park Rd. Schaumburg, IL 60193 U.S.A. 1-847-534-6000 1-800-247-3368 www.bisco.com

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