BISCO

CE 0459

VELITEFLO™ LV

Light-Cured

Radiopaque Low Viscosity Composite

Instructions for Use

BISCO, Inc.
1100 W. Irving Park Rd.
Schaumburg, IL 60193
U.S.A.
1-847-534-6000
1-800-247-3368

IN-1947R
Rev. 5/22

en	VELITEFLO™ LV Radiopaque Low Viscosity Composite
es	VELITEFLO™ LV Resina compuesta radiopaca de baja viscosidad
fr	VELITEFLO™ LV Composité radio-opaque de faible viscosité
de	VELITEFLO™ LV Strahlennundschüssiges Komposit mit niedriger Viskosität
it	VELITEFLO™ LV Composito radiopaco a bassa viscosità
sv	VELITEFLO™ LV Radiopak, lågvisköst komposit
da	VELITEFLO™ LV Røntgenpositivt komposit af lav viskositet
fi	VELITEFLO™ LV Röntgenpositiivinen pienviskoinen -komposiitti
nl	VELITEFLO™ LV Radiopake laagviskeuze composit
no	VELITEFLO™ LV radiopak komposit med lav viskositet
pt	VELITEFLO™ LV Composito radiopaco de baixa viscosidade
el	VELITEFLO™ LV Ακτινοεπαρκές σύνθετο ρητίνης χαμηλού βύσους
ru	VELITEFLO™ LV Аκτινοεπαρκές σύνθετο ρητίνης χαμηλού βύσους
bg	VELITEFLO™ LV Рентгеноконтрастен композитен материал с ниска вискозитет

VELITEFLO LV is a light-cured, highly radiopaque microhybrid composite with a low modulus of elasticity.
Indications For Use:

- Cavity liner (Class I or Class II)
- Pit and fissure sealant
- Preventative Resin Restorations
- Blockout undercuts

Warnings:

- If accidentally splashed into the eye, flush with copious amounts of water and seek medical attention immediately. In case of contact with other tissues, rinse immediately with plenty of water for several minutes.

Intended user: Licensed dentist.
Caution:

- Product may contain items that are designed for one time use. Dispose of used or contaminated tips. Do not clean, disinfect or reuse.
- Use of commonly available hygienic protective covering to avoid any contamination of syringes during treatment is recommended. If syringe becomes contaminated, discard. Do not clean or disinfect.

Precautions:

- Avoid contact with the skin; unpolymerized (meth)acrylate resins may cause skin sensitization in susceptible persons. In case of contact, wash with adhesive soap and water.
- When using dental adhesives, saliva contamination will seriously compromise dentin bonding.
- Refer to individual component labels for specific expiration dates.
- Protective eyewear should be worn by the patient and dental staff.
- See individual component labels for specific expiration dates.
- Safe data sheet available on request.
- Safety data sheet available at www.bisco.com.
- Instructions for Use available on request.
- Instructions for Use available at www.bisco.com.

Helpful Hints:

- The restoration may appear smooth and polished immediately after placement, however an oxygen-inhibited layer is present and must be removed to prevent staining.

INSTRUCTIONS FOR USE

- Prior to isolation, select appropriate shade of **VELITEFLO LV**.
- Remove the cap from the selected syringe and attach a needle tip securely.
- Be sure the cavity is adequately cleaned in accordance with the Instructions for Use of the bonding system to be used.
- Apply an adhesive according to the manufacturer's instructions.
- Using the previously selected shade, place 1-2mm increments of composite into the cavity preparation.
- Light cure each increment for the 20 seconds. Continue to build incrementally as desired.

HYGIENE: Re-cap syringes with the Luer lock cap. Use of commonly available hygienic protective covering to avoid any contamination of syringes during treatment is recommended.
DISPOSAL: Refer to community provisions relating to waste. In their absence, refer to national or regional provisions relating to waste.
STORAGE: Store at room temperature (20°/68°F - 25°/77°F).
WARRANTY: BISCO, Inc. recognizes its responsibility to replace products if proven to be defective. BISCO, Inc. does not accept liability for any damages or loss, either direct or consequential, stemming from the use of or inability to use the products as described. Before using the product, the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.
NOTICE: Report any serious incident to BISCO and the competent authority of the state in which the user and/or patient is established.
Summary of Safety and Clinical Performance can be found at <https://ec.europa.eu/tools/edumad> or available upon request.
To Reorder Accessories:
X-80608N – 50 Blue Disposable Tips (22 Gauge)
* VELITEFLO is a trademark of BISCO, Inc.

en	VELITEFLO™ LV Radiopaque Low Viscosity Composite
es	VELITEFLO™ LV Resina compuesta radiopaca de baja viscosidad
fr	VELITEFLO™ LV Composité radio-opaque de faible viscosité
de	VELITEFLO™ LV Strahlennundschüssiges Komposit mit niedriger Viskosität
it	VELITEFLO™ LV Composito radiopaco a bassa viscosità
sv	VELITEFLO™ LV Radiopak, lågvisköst komposit
da	VELITEFLO™ LV Røntgenpositivt komposit af lav viskositet
fi	VELITEFLO™ LV Röntgenpositiivinen pienviskoinen -komposiitti
nl	VELITEFLO™ LV Radiopake laagviskeuze composit
no	VELITEFLO™ LV radiopak komposit med lav viskositet
pt	VELITEFLO™ LV Composito radiopaco de baixa viscosidade
el	VELITEFLO™ LV Ακτινοεπαρκές σύνθετο ρητίνης χαμηλού βύσους
ru	VELITEFLO™ LV Аκτινοεπαρκές σύνθετο ρητίνης χαμηλού βύσους
bg	VELITEFLO™ LV Рентгеноконтрастен композитен материал с ниска вискозитет

VELITEFLO LV is a light-cured, highly radiopaque microhybrid composite with a low modulus of elasticity.
Indications For Use:

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