

TheraCal PT™

Bisco Dental
www.bisco.com



Clinician: Dr. Ona Erdt

Problem:

Wiggly, salivating pediatric patients having pulpal treatment must sit for several minutes with medicaments in the tooth. In addition, other pulpotomy materials may require hand mixing a powder and a liquid and often have a long set time.



Solution:

TheraCal PT™ is a resin-modified, dual-cured, biocompatible material designed for pulpal treatment that is light cured for 10 seconds. This allows for an immediate or same-day restoration. Once hemostasis is achieved, **TheraCal PT** can be conveniently placed from an automix syringe. In addition to saving time, **TheraCal PT** is highly radiopaque, releases calcium and has an alkaline pH.

