

# 5 REASONS TO TRY TheraCem<sup>®</sup> Ca

Dual-Cured, Self-Adhesive Calcium Releasing Cement

Rx Only

## 1 Contains Calcium

Continuous calcium ion release.\*

## 2 Contains MDP

Strong bond to zirconia, metal and alumina substrates without primer.

## 3 Easy Clean-up

Tack curing for 2-3 seconds facilitates the removal of excess cement.

## 4 Alkaline pH

Transitions from acidic to an alkaline pH in minutes.\*

## 5 Radiopaque

Easily identified on radiographs for quick and effective diagnosis.



\*Data on file. BISCO, Inc.

MC-50231TCI

For more information visit [www.bisco.com](http://www.bisco.com)

