

UNIVERSALLY  
**STRONG.**  
UNIVERSALLY  
**SIMPLE.**





SCAN TO  
LEARN MORE!



## CASE STUDY

Courtesy of Adamo Notarantonio, DDS, AAACD, FICOI



## ORDERING INFORMATION

**DUO-LINK UNIVERSAL SYSTEM KIT**..... A-19620K  
 1 Dual-Syringe DUO-LINK UNIVERSAL Shade: UNIVERSAL (8g), 1 Bottle ALL-BOND UNIVERSAL (6ml),  
 1 Bottle Z-PRIME PLUS (2ml), 1 Bottle Porcelain Primer (3ml), Accessories, Instructions

### 1 DUAL-SYRINGE PACKAGE

1 Dual-syringe (8g), 15 Mixing Tips and Intraoral Tips, Instructions

**UNIVERSAL** ..... A-19030P  
**MILKY WHITE** ..... A-197MWP

**ALL-BOND UNIVERSAL (6ml Bottle)** ..... B-7202P

**Z-PRIME PLUS (2ml Bottle)** ..... B-6002P

**ZIRCLEAN® (5g Syringe)** ..... B-7351P

**PORCELAIN PRIMER (10ml Bottle)** ..... B-2221P

Call us! We're here to help:

1-800-247-3368

www.bisco.com



# DUO-LINK UNIVERSAL™



# DUO-LINK UNIVERSAL™

The adhesive resin cement is specially formulated for cementation of ALL† indirect restorations. It is available in two esthetic shades, Universal and Milky White. Duo-Link Universal is intended for use with adhesives designed for compatibility with all dental materials, including all BISCO adhesives.



## THE PERFECT FIT FOR YOUR OFFICE



**Extremely easy clean-up**  
Formulated to allow for quick and easy removal of excess cement



**Highly radiopaque**  
Easily identified on radiographs for quick and effective diagnosis



**CAD/CAM restorations**  
Ideal for all chairside and lab-fabricated restorations




**Universal for all cementation procedures†**  
Indicated for all crowns, bridges (fixed prosthesis), inlays, onlays, and posts/dowels, fabricated from metal, composite, porcelain, glass ceramics, zirconia, and alumina.



**Easy to use**  
Auto-mix, dual-syringe provides a consistent mix for immediate delivery



**High degree of conversion**  
High degree of conversion in both light- and self-cured modes ensures a strong, long lasting restoration



**Two shades**  
Universal shade is a translucent A2 and the Milky White shade is more opaque and can be used to mask darker tooth preps

† It is recommended to use BISCO's CHOICE™ 2 for veneer cementation.

## WHAT'S IN THE KIT?

The Duo-Link Universal kit is an adhesive cementation system. The combination of award-winning, innovative products provides high retention and durability of indirect substrates in a simple, universal, and easy-to-use kit.

### DUO-LINK UNIVERSAL

This adhesive resin cement specially formulated for cementation of ALL† indirect restorations

### ALL-BOND UNIVERSAL

Light cure, MDP containing, truly universal, single bottle adhesive

### PORCELAIN PRIMER

Silane coupling agent for all glass ceramics surfaces


### Z-PRIME™ PLUS

Single component and MDP containing primer for zirconia and metal surfaces




## SEE WHAT THE DATA SAYS


### SBS Zirconia (Sandblasted)



### SBS e.max (Etched)




### SBS Dentin (Total Etch)




\*All cements were used in light-cured mode. Trademarks are property of their respective owners.

# DUO-LINK UNIVERSAL PREP AND CEMENTATION STEP BY STEP

## 1 RESTORATION TREATMENT\*\*




## 2 ETCHING



## 3 ADHESIVE APPLICATION



## 4 CEMENTATION



\*\* Refer to complete instructions for use for detailed technique and instructions.